



SKORR | GLOW

RECOMMENDED TREATMENT PROTOCOL

Choose the appropriate mode based on your treatment goals and desired results. Always follow an upward and outward motion when using the device. Begin treatment with 3 gliding passes and complete treatment with 1-3 passes in stamping mode.

RECOMMENDED FREQUENCY

Use each mode 1-2 times daily for optimal results.

PRECAUTIONS

- Avoid use immediately after medical procedures or surgery.
- Do not use on open wounds, infected, or inflamed skin.
- Consult a physician prior to use, if needed.
- For first-time users, start with the lowest intensity setting and adjust gradually based on comfort.



LEARN MORE

www.skorrglow.com

User Guide

SKORR
DR. ANETTA RESZKO

HOW TO USE **SKORR** | GLOW



BOOST Mode

Technology: Electroporation
+ Yellow Multiwave LED

Benefits:

- Maximize absorption of skincare ingredients and help deliver intense hydration
- Promote lymphatic drainage to reduce puffiness
- Soothe and calm skin to minimize inflammation



LIFT Mode

Technology: Sculpting & Toning Microcurrent
+ Green Multiwave LED

Benefits:

- Firm and revitalize skin for a youthful appearance
- Sculpt facial contours to achieve a radiant, lifted and snatched look
- Smooth fine lines and wrinkles
- Reduce puffiness and inflammation by promoting lymphatic drainage



BUILD Mode

Technology: Rejuvenating Radiofrequency
+ Sonic Vibrations + Red Multiwave LED

Benefits:

- Boost collagen production for restored skin elasticity and firmness
- Visibly reduce fine lines and wrinkles
- Minimize the appearance of pores
- Improve circulation for enhanced oxygen and nutrient delivery
- Support skin detoxification for a clearer, healthier complexion



EYE & LIP Mode

Technology: Rejuvenating Radiofrequency
+ Sonic Vibrations

Benefits:

- Hydrate and plump lips
- Depuff and rejuvenate eyes
- Diminish fine lines, wrinkles, and sagging
- Brighten and even out discoloration for a more luminous appearance